XPUR

M-70-15-1316 / PSD-0000058 (04) / JUL - 2022

Instructions for use
ENGLISH

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Promedon

DESCRIPTION

The structure of the xPur dilation balloons can be divided into:

Body of the balloon, located in the proximal end, has two radiopaque marker bands that show the nominal length of the balloon. The balloon has a stable shape and allows high operational pressure. (Nominal pressure shown on label). The working pressure (nominal pressure) exceed the recommended limit, due to the risk of balloon rupture.

Luer-lock connectors, located in the distal end, allow external link of the balloon using the catheter. The catheter has two lumens which can be accessed from the luer-lock connectors. The central lumen runs through the device linking both ends and is used for concentric access of the 0.035" guide wire. The other lumen connects to the interior of the balloon and it is used to increase its internal pressure in order to reach the required size.

xPur ureteral dilation balloons have a protective cover which must be peeled off before its use, as well as a metallic guide protecting the central lumen of the catheter, which must be removed in order to place the guide wire. The dimensions of the Balloon's body, as well as the catheter, are shown on the product's label.

CONTENTS

xPur comprises:

1 (one) ureteral dilation balloon. (Size indicated on label)

The product is protected by a plastic tube, inside a protective envelope.

INDICATIONS

PROMEDON's xPur ureteral dilation balloon has been developed for the dilation of narrowed and stenotic ureters.

WARNINGS

Over dilation of the affected section must be avoided. The vessel's diameter should not exceed the pre-stricture segment's diameter.

PRECAUTIONS

Single-use sterile product (sterilized with ethylene oxide); must not be reused or resterilized.

Do not use the kit, if the package is open or damaged.

It is advisable to store the kit in its original packaging at: Room TEMPERATURE.

USE BEFORE THE EXPIRY DATE SPECIFIED ON THE PACKAGING.

Only a physician with experience in the application of the kit is qualified to use it.

CONTRAINDICATIONS

- Stenosis that cannot be crossed by the catheter.
- Patients whose injury prevents full inflation of the balloon.

COMPLICATIONS

Ureteral dilation procedures should only be performed by physicians familiar with the inherent potential complications. The entire process must be carried out under radiographic control. Complications may occur during or after the procedure, and may include the following:

- Dissection.
- Perforation of the ureter.
- Hematomas or hemorrhages.
- Allergy to the contrast substance.
- Infections or sepsis.

PREPARATION

- Make sure the product's packaging is not damaged. Do not use the product if the packaging is damaged.
- Before the procedure, remove the product from its packaging in a suitable area, taking the necessary precautions for handling sterile material.
- Verify that the product is not damaged. Damaged products must not be used.
- Carefully remove the balloon's protection from the proximal end of the catheter. Avoid pulling from the dilation balloon's ends.

Deflating the dilation balloon

- Connect a 2 or 3 way stopcock (not included) to the dilation balloon's luer-lock inflation connection.
- Connect a syringe and open the 2 or 3 way stopcock. Make sure that the syringe volume is adequate for the type of dilation balloon being used.
- Draw back on the syringe, applying negative pressure to the balloon.
- Flush the central lumen (through where the guide is inserted) with saline solution.

- Carefully take the dilation balloon catheter, using the guide wire previously placed, through the stenotic ureter's lumen.
- Place the dilation balloon under radiographic control in the dilation area. The marker bands are used as a reference guide for correct location of the catheter, and they define the nominal length of dilation.
- The guide wire must be advanced beyond the distal end of the catheter.

Inflating the dilation balloon

- Balloon inflation should be carried out using a syringe or a manometer.
- Connect a syringe (or manometer) to the 2 or 3 way stopcock.
- Inject the liquid to inflate the dilation balloon, slowly and evenly. Maintain the pressure for a few seconds. Do not exceed the nominal rupture pressure.

Deflating the dilation balloon

- To deflate the dilation balloon, draw a vacuum with a syringe. In the event that a liquid is used, slowly draw the syringe's piston back, and maintain the low pressure for a while to allow the complete evacuation of the liquid.
- Slowly remove the deflated dilation balloon, twisting slightly, maintaining the guide wire in place.

Indications for proper use of the dilation balloon.

- It is advised to use a manometer to maintain the inflation pressure both constant and visible.
- If there is resistance on removing the dilation balloon, the entire device should be carefully removed as a unit to avoid damaging the ureter.
- Should loss of pressure occur when inflating the dilation balloon, caused by rupture, the dilation balloon should be carefully removed using the guide wire.

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Fabricado por: PROMEDON S.A
PROMEDON S.A.