




introducing
the **OSV II**[®] and
OSV II Low Flow[™]
lumbar valve systems

The only shunts that **actively manage** patients' changing drainage needs **without** the need for time consuming adjustments.




introducing the **OSV II**[®] and **OSV II Low Flow**[™] lumbar valve systems

- **No time-consuming adjustments**
- **MR Safe**—composed of all non-ferromagnetic material
- **Simple and quick assembly reduces procedure time**
- **Ensures a close balance between CSF production and drainage rate**



The OSV II[®] achieves these feats with its patented flow regulating technology. It is able to transition from a valve with a low opening pressure to one with a higher pressure resistance by reducing the size of the valve opening as pressure increases. During this scenario, the valve works to maintain a flow rate of approximately 20 ml/hr (similar to the natural rate of production in adults).



Other lumbar systems simply open and close based on pressure differential; however, **these valves fail to account for CSF flow rate**. As a result, whenever there is a pressure increase derived from hydrostatic or vasogenic causes; such as REM sleep, standing up, or increased physical activity; a typical lumbar shunt system


A Physiological Approach to Lumbar Shunting that Reduces the Likelihood of Overdrainage

will allow fluid to drain at excessive rates until pressure is relieved. It is this rapid drainage that can increase adverse effects such as over-drainage and proximal catheter obstruction.


The OSV II system was designed to mimic normal physiological behavior by automatically adjusting CSF flow levels in response to differential pressure and ensuring that CSF flow does not reach excessive rates. Maintaining a close balance between CSF production and drainage rate allows the ventricles to reduce gradually.

This innovative technology is now available for implantation in the lumbar subarachnoid space. The complete procedure ready kit requires one simple connection, therefore reducing assembly time in the OR suite.

In essence the OSV II system is a unique, easy to use system that aims to mimic normal physiological behavior by maintaining steady CSF flow rates. Instead of patients having to constantly adjust to their shunt system, the system constantly adjusts to them.




By adjusting the size of the valve passageway, the OSV II system works to maintain a flow rate of approximately 20 ml/hr (consistent with the natural rate of production in adults).



For patients that require a drainage rate closer to 10 ml/hr*, Integra offers the OSV II Low Flow Lumbar Valve System.

*CSF production is reduced as part of healthy aging. This requires that some patients maintain lower flow rates.



U.S. Ordering Information:

- 909722:** OSV II Lumbar Valve System with Antechamber
- 909723:** OSV II Lumbar Valve System without Antechamber
- 909522:** OSV II Low Flow Lumbar Valve System with Antechamber
- 909523:** OSV II Low Flow Lumbar Valve System without Antechamber

International (outside of the United States) Ordering Information:

- 909710:** OSV II Lumbar Valve System with Antechamber
- 909711:** OSV II Lumbar Valve System without Antechamber
- 909510:** OSV II Low Flow Lumbar Valve System with Antechamber
- 909511:** OSV II Low Flow Lumbar Valve System without Antechamber

Kit Contains:



- 80 cm, F5 closed tip lumbar catheter
- Valve unit with integral open ended 110 cm F7 peritoneal catheter
- 14G, 9 cm Touhy Needle with snap on wings
- PTFE coated guidewire in dispenser
- Polypropylene stepdown connector (F8 to F5)
- Polypropylene straight connector
- Suturable tubing clamp
- Luer lock connector



OPERATE WITH CONFIDENCE™



INDICATIONS:

The Integra NPH™ Low Flow Valve is an implantable system used in the treatment of patients with hydrocephalus, to shunt CSF from the ventricles to the peritoneal cavity or other appropriate drainage site such as the heart's right atrium.

The OSV II® System is an implantable system used in the treatment of patients with hydrocephalus, to shunt CSF from the ventricles to the peritoneal cavity or other appropriate drainage site such as the heart's right atrium.

Integra LifeSciences Corporation
311 Enterprise Drive
Plainsboro, NJ 08536
www.Integra-LS.com

USA and Canada:
800.654.2873
609.275.0500 (Outside USA)
609.275.5363 (Fax)

England:
+44 (0) 1264 345 700
+44 (0) 1264 332 113 (Fax)

France:
+33 (0) 493 95 56 00
+33 (0) 493 95 56 60 (Fax)

Germany:
0800 10 10 755
06995 775 477 (Fax)

Camino, CUSA, DuraGen Plus, Integra NeuroSciences, OSV II, and Radionics are registered trademarks of Integra LifeSciences Corporation or its subsidiaries. Auragen, Integra NPH, Operate with Confidence, Ruggles, and the Integra wave logo are trademarks of Integra LifeSciences Corporation or its subsidiaries. MAYFIELD is a registered trademark of SM USA Inc. and is used by Integra LifeSciences Corporation under license.
©2010 Integra LifeSciences Corporation. All rights reserved. Printed in the USA. NS3800-02/10

